




# NeoChild<sup>®</sup>

(888) 887-NICU (6428)  
www.neochild.com

## NeoChild Umbilical Catheterization Tray



### NeoChild Forceps 2.75"

(Order #: NC20001,  
10 Forceps per Case)

### Tray Contents (Order #NC20000T: 5 Trays per Case)

- |                       |                     |                                 |
|-----------------------|---------------------|---------------------------------|
| (1) Adson Forceps     | (1) Scissors        | (1) Filter Needle               |
| (1) Iris Forcep 4"    | (1) 3-Way Stopcock  | (1) 3-Pack Iodophor Swab        |
| (1) Iris Forcep 2"    | (1) 1cc Syringe     | (1) CSR Wrap                    |
| (1) Needle Holder     | (1) 3cc LL Syringe  | (2) 18" x 16" Poly-Lined Drape  |
| (1) Scalpel #11 Blade | (1) Tape Measure    | (1) 12" x 18" Fenestrated Drape |
| (1) Hemostat Straight | (1) 5" Probe        | (4) 4" x 4" Gauze Pads          |
| (1) 4-0 Silk Suture   | (1) Sodium Chloride | (4) 2" x 2" Gauze Pads          |
| (1) Hemostat Curved   | (1) 20ga Needle 1½" | (2) Point-Lok Sharps Container  |



NeoChild<sup>®</sup>

(888) 887-NICU (6428)  
www.neochild.com