



# NeoChild®

Creator of the  
Safe Child System™

(888) 887-NICU (6428)  
www.neochild.com

## Enteral Feeding Administration Sets

### EXTENSION SETS

\*\*\*All with orange striping and tagged "For Enteral Feedings Only"

**IS18ENO** 18-inch tubing, orange female enteral and male enteral adapter (non-IV).

**IS36ENO** 36-inch tubing, orange female enteral and male enteral adapter (non-IV).

**IS60ENO** 60-inch tubing, orange female enteral and male enteral adapter (non-IV).

**IS60ENYO-C** 60-inch tubing with 'Y', 2 orange female enteral adapters and 1 male enteral adapter (non-IV).

**IS36EN** 36-inch tubing and orange male enteral only adapter.

**IS60EN** 60-inch tubing and orange male enteral only adapter.

**IS60ENY-C** 60-inch tubing with 'Y' and an orange male enteral only adapter.

### SILICONE FEEDING TUBES – ENTERAL ONLY

	<u>40 cm (16")</u>	<u>60 cm (24")</u>	<u>90 cm (36")</u> *
5.0 FR	NC1S50SE	NC1S50ME	NC1S50LE*
6.5 FR	NC1S65SE	NC1S65ME	NC1S65LE*
8.0 FR	NC1S80SE	NC1S80ME	NC1S80LE*

\*Special order feeding tubes, contact your representative for assistance

### POLYURETHANE FEEDING TUBES – ENTERAL ONLY

	<u>40 cm (16")</u>	<u>60 cm (24")</u>	<u>90 cm (36")</u>
4.0 FR	NC1P40SE	NC1P40ME	
5.0 FR	NC1P50SE	NC1P50ME	NC1P50LE
6.5 FR	NC1P65SE	NC1P65ME	NC1P65LE
8.0 FR	NC1P80SE	NC1PB0ME	NC1P80LE

### PVC FEEDING TUBES – ENTERAL ONLY

	<u>40 cm (16")</u>	<u>90 cm (36")</u>
5.0 FR	NC1V50SE	NC1V50LE
6.0 FR	NC1V60SE	NC1V60LE
8.0 FR	NC1V80SE	NC1V80LE

### ENTERAL SYRINGES – ORANGE-TIPPED ENTERAL ONLY

<u>Size</u>	<u>Order No.</u>	<u>No. / Case</u>
1 mL	NC600100	100
3 mL	NC600300	100
5 mL	NC600500	100
10 mL	NC601000	100
20 mL	NC602000	100
30 mL	NC603000	50
60 mL	NC606000	25



NeoChild®

(888) 887-NICU (6428)  
www.neochild.com