




NeoChild[®]

Safe Child System[™]
US Patent 8,292,875

(888) 887-NICU (6428)
www.neochild.com


The Safe Child System[™]




NeoChild's innovative SafeChild System[®] eliminates enteral tubing misconnections and provides you with a safe and high-quality feeding system.



Enteral-only BD syringes



Enteral extension sets



Polyurethane and Silicone feeding tubes



NeoChild[®]

(888) 887-NICU (6428)
www.neochild.com



NeoChild[®]

Safe Child System[™]
US Patent 8,292,875

(888) 887-NICU (6428)
www.neochild.com

Enteral Feeding Administration Sets

EXTENSION SETS

***All with orange striping and tagged "For Enteral Feedings Only"

IS18ENO 18-inch tubing, orange female enteral and male enteral adapter (non-IV).

IS36ENO 36-inch tubing, orange female enteral and male enteral adapter (non-IV).

IS60ENO 60-inch tubing, orange female enteral and male enteral adapter (non-IV).

SILICONE FEEDING TUBES—ENTERAL ONLY

	<u>40 cm (16")</u>	<u>60 cm (24")</u>	<u>90 cm (36")*</u>
5.0 FR	NC1S50SE	NC1S50ME	NC1S50LE*
6.5 FR	NC1S65SE	NC1S65ME	NC1S65LE*
8.0 FR	NC1S80SE	NC1S80ME	NC1S80LE*

*Special order feeding tubes, contact your representative for assistance

POLYURETHANE FEEDING TUBES—ENTERAL ONLY

	<u>40 cm (16")</u>	<u>60 cm (24")</u>	<u>90 cm (36")</u>
4.0 FR	NC1P40SE	NC1P40ME	
5.0 FR	NC1P50SE	NC1P50ME	NC1P50LE
6.5 FR	NC1P65SE	NC1P65ME	NC1P65LE
8.0 FR	NC1P80SE	NC1PB0ME	NC1P80LE

BD ENTERAL SYRINGES—ENTERAL ONLY

<u>Size</u>	<u>Order No.</u>
1 mL	305851
3 mL	305853
5 mL	305855
10 mL	305857
20 mL	305859
30 mL	305861
60 mL	305863



NeoChild[®]

(888) 887-NICU (6428)
www.neochild.com