



NeoChild[®]

Creator of the
Safe Child System[™]

(888) 887-NICU (6428)
www.neochild.com

Urinary Drainage Closed Collection Kits



Tray Contents

- (1) 150mL Collection Burette
- (1) 5-Ply Poly Drape Towel
- (1) 3-Pack Iodophor Swabs
- (1) NeoCath Drape
- (4) 2" x 2" Gauze
- (1) Catheter Lubricant



Order Information (5 Trays per Case):

Size	No Catheter	3.5 FR	5.0 FR	6.5 FR	8.0 FR
Order #	NC30000T	NC30350T	NC30500T	NC30650T	NC30800T



NeoChild[®]

(888) 887-NICU (6428)
www.neochild.com