





Infusion Devices™

(888) 887-NICU (6428)
www.neochild.com

Multiple Access IV Extension Sets



Order #: IS6003 (50 per Case)
Microbore Extension Set.
60", 0.3 mL priming volume.
Male/Female luer lock and removable slide clamp.



Order #: IS6010 (50 per Case)
Minibore Extension Set.
60", 1.0 mL priming volume.
Male/Female luer lock and pinch clamp.



Order #: IS3618 (50 per Case)
Microbore Extension Set.
36", 0.18 mL priming volume.
Male/Female luer lock and pinch clamp.



Order #: IS6010SC (50 per Case)
Minibore Extension Set.
60", 1.0 mL priming volume.
Male/Female luer lock, slide clamp and stopcock.



Order #: IS7220 (50 per Case)
Minibore Extension Set.
72", 2.0 mL priming volume.
Male/Female luer lock and slide clamp.



Order #: IS6018 (50 per Case)
Minibore Extension Set.
60", 1.8 mL priming volume.
Male luer lock, needleless valve and pinch clamp.

**For additional sets, please contact
info@neochild.com**



NeoChild®
(888) 887-NICU (6428)
www.neochild.com