



NeoChild[®]

PRODUCT CATALOG



NeoChild[®]
(888) 887-NICU (6428)
www.neochild.com




TABLE OF CONTENTS

Enteral Only Syringes	2
Enteral Only Extension Sets	2
Enteral Only Feeding Tubes (Silicone)	3
Enteral Only Feeding Tubes (Polyurethane)	3
Umbilical Catheterization Trays	4
Stainless Forceps	4
PICC Insertion Trays	4
Urinary Drainage Catheters (Silicone)	5
Urinary Drainage Trays	5
Lumbar Puncture Needles	5
Lumbar Puncture Trays	5
IV Extension Sets	6–9
Hydrogel Probe Covers	9



The Safe Child System

Eliminate misconnections with the safety of enteral only connections with feeding tubes, extension sets, and syringes.




■ Enteral Syringes

Enteral only syringes (legacy, non-IV).

Size	Order #
1 mL Syringe	NC600100
3 mL Syringe	NC600300
5 mL Syringe	NC600500
10mL / 12 mL Syringe	NC601000 / NC601200
20 mL Syringe	NC602000
30 / 35 mL Syringe	NC603000 / NC603500
60 mL Syringe	NC606000

Oral tip cap included




■ Enteral Only Extension Sets

Enteral only tubing with orange female and male connectors (legacy, non-IV).

Size	Order #	Case Size
1 mL (18")	IS18ENO *	50 per Case
1.5 mL (36")	IS36ENO *	50 per Case
2 mL (60")	IS60ENO	50 per Case

*Bifurcated options available upon request; * Indicates Special Order*

**** 5" Enteral Only Sampling Straws Also Available (NCSTRAW)**



▪ Silicone Enteral Only Feeding Tubes

Enteral Only Feeding Tube with Orange Radiopaque Stripe and Enteral Female Connector.

Size	Order #	Case Size
5.0 FR (40cm/Short)	NC1S50SE	10 per Case
5.0 FR (60cm/Medium)	NC1S50ME	10 per Case
5.0 FR (90cm/Long)	NC1S50LE	10 per Case
6.5 FR (40cm/Short)	NC1S65SE	10 per Case
6.5 FR (60cm/Medium)	NC1S65ME	10 per Case
6.5 FR (90cm/Long)	NC1S65LE	10 per Case
8.0 FR (40cm/Short)	NC1S80SE	10 per Case
8.0 FR (60cm/Medium)	NC1S80ME	10 per Case
8.0 FR (90cm/Long)	NC1S80LE	10 per Case

Special Order Item, please contact representative for details

▪ Polyurethane Enteral Only Feeding Tubes

Enteral Only Feeding Tube with Orange Radiopaque Stripe and Enteral Female Connector.



Size	Order #	Case Size
4.0 FR (40cm/Short)	NC1P40SE	10 per Case
4.0 FR (60cm/Medium)	NC1P40ME*	10 per Case
5.0 FR (40cm/Short)	NC1P50SE	10 per Case
5.0 FR (60cm/Medium)	NC1P50ME	10 per Case
5.0 FR (90cm/Long)	NC1P50LE	10 per Case
6.5 FR (40cm/Short)	NC1P65SE	10 per Case
6.5 FR (60cm/Medium)	NC1P65ME	10 per Case
6.5 FR (90cm/Long)	NC1P65LE	10 per Case
8.0 FR (40cm/Short)	NC1P80SE	10 per Case
8.0 FR (60cm/Medium)	NC1P80ME	10 per Case
8.0 FR (90cm/Long)	NC1P80LE	10 per Case

** Indicates Special Order*



NeoChild®

info@neochild.com

(888) 887-NICU (6428)

www.neochild.com

ISO-80369 Compatible System

FEEDING TUBES	40cm (16")	60cm (24")	90cm (36")	SYRINGES	With Oral Tip Cap
5.0FR	NC1P50SE-ENF	NC1P50ME-ENF	NC1P50LE-ENF	1ml	NC600100-ENF
6.5FR	NC1P65SE-ENF	NC1P65ME-ENF	NC1P65LE-ENF	3ml	NC600300-ENF
8.0FR	NC1P80SE-ENF	NC1P80ME-ENF	NC1P80LE-ENF	5ml	NC600500-ENF
				10ml	NC601000-ENF
				20ml	NC602000-ENF
				30ml	NC603000-ENF
				60ml	NC606000-ENF

EXTENSION SETS	2.0ml Approximate Priming Volume
60"	IS60ENO-ENF
60" with Medication Port	IS60ENYO-ENF



NeoChild's ISO 80369-3 compliant system offers on-shore manufacturing to solidify supply chains and offer competitive price points.



ISO-80369 Syringes



ISO 80369-3 Extension Sets



ISO 80369-3 Feeding Tubes

** Additional Accessories Available **




NeoChild®



▪ NeoChild Umbilical Catheterization Tray

Complete tray (no catheter) for umbilical vessel catheter placement.

Tray Contents— Order #: NC20000T (10 Per Case)



- (1) Adson Forcep, 4.75"
- (1) Point-Lok™ Sharps Container
- (1) Needle Holder, 5"
- (1) Safety Scalpel
- (1) 4-Way Stopcock
- (1) Straight Hemostat, 5.5"
- (2) Curved Forceps, 5"
- (1) Scissors, 5.5"
- (1) 1cc Syringe TB L/S
- (1) 5cc Syringe L/L
- (2) 3cc Syringe L/L
- (1) Tape Measure
- (1) Probe, 5"
- (1) 3-Pack PVP Swabsticks
- (1) 20ga x 1" Needle
- (1) 18ga x 1.5" Filter Needle
- (1) Curved Forceps, 3.5"
- (1) CSR Wrap
- (1) Umbilical Tie,
- (2) 18" x 26" Polylined Drape Towel
- (1) Fenestrated Drape
- (1) 4-0 Silk Suture
- (4) Gauze Pads, 4" x 4"
- (4) Gauze Pads, 2" x 2"

▪ Stainless Forceps 3.5"

Stainless steel forceps for umbilical vessel catheter placement.



Size	Order #	Case Size
3.5"	NC200001	10 per Case

▪ NeoChild PICC Insertion Tray

Complete tray (no catheter) for peripherally inserted central catheter placement.

Tray Contents— Order #: NC50000T (20 Per Case)



- (2) Tape Measure
- (4) 2" x 2" 12-Ply Gauze Pads
- (1) PVP Swabsticks (Triple)
- (1) Saline Wipe
- (3) Steri Strip
- (1) Fenestrated Drape
- (1) CSR Wrap
- (1) Adson Forceps, 4.75"
- (1) Curved Forceps, 3.5"
- (1) Transpore Tape
- (1) Tegaderm
- (1) Poly-Lined Drape
- (1) Vial Spike
- (1) Silicone Tourniquet
- (1) Scissors, 5.5"
- (1) 10cc LL Syringe
- (6) 4" x 4" 12-Ply Gauze Pads
- (1) Skin Prep Pad
- (3) Labels



* Indicates Special Order

NeoChild Silicone Urinary Drainage Catheters Radiopaque Stripe.

Size	Order #	Case Size
3.5 FR (40cm)	NC303500 *	10 per Case
5.0 FR (40cm)	NC305000	10 per Case
6.5 FR (40cm)	NC306500	10 per Case
8.0 FR (40cm)	NC308000	10 per Case



NeoChild Urinary Drainage Closed Collection Kits Complete kit with or without silicone catheter.

Size	Order #	Case Size	Tray Contents
No Catheter	NC30000T *	10 per Case	(1) 100 ml Collection Burette (1) Antiseptic Swabsticks (1) Catheter Lubricant (1) Polylined Drape (1) Fenestrated Drape (4) Gauze Pads

Lumbar Puncture Needles




Size	Order #	Case Size
22ga x 1"	NC401000	10 per Case
22ga x 1.5"	NC401500	10 per Case

NeoChild Lumbar Puncture Kits




Size	Order #	Case Size	Tray Contents
22ga x 1" Tray	NC40100T	5 per Case	(1) 22ga x 1" or 1.5" Spinal Needle
22ga x 1.5" Tray	NC40150T	5 per Case	(4) Numbered Specimen Tube with Screw Cap (1) CSR Wrap (1) 4-Ply Poly Drape Towel (1) Point-Lok Sharps Container (1) 3-Pack PVP Swabsticks (1) Spot Adhesive Bandage (1) Neocath Drape (4) 2" x 2" Gauze



- **Order #: IS0702* (50 per Case)**
Microbore Extension Set.
7", 0.2 mL priming volume.
Male/Female luer lock and slide clamp.



- **Order #: IS3618 (50 per Case)**
Microbore Extension Set.
36", 0.18 mL priming volume. Male/Female luer lock and pinch clamp.



- **Order #: IS6003 (50 per Case)**
Microbore Extension Set.
60", 0.3 mL priming volume. Male/Female luer lock and removable slide clamp.




- **Order #: IS6018 (50 per Case)**
Minibore Extension Set.
60", 1.8 mL priming volume. Male luer lock, needleless valve and pinch clamp.



- **Order #: IS6010 (50 per Case)**
Minibore Extension Set.
60", 1.0 mL priming volume. Male/Female luer lock and slide clamp.

**Indicates special order, please contact representative for details*



▪ **Order #: IS6010SC (50 per Case)**

Minibore Extension Set.
60", 1.0 mL priming volume. Male/Female luer lock, slide clamp and stopcock.




▪ **Order #: IS6018DT (50 per Case)**

Minibore Dark Tubing Extension Set for light sensitive medications.
60", 1.8 mL priming volume. Male/Female luer lock and pinch clamp.



▪ **Order #: IS72V* (50 per Case)**

Minibore Extension Set with Anti-siphon Valve.
72", 2.0 mL priming volume. Male/Female luer lock and removable slide clamp.



▪ **Order #: IS7220 (50 per Case)**


Minibore Extension Set.
72", 2.0 mL priming volume. Male/Female luer lock and slide clamp.

**Indicates special order, please contact representative for details*




▪ **Order #: IS3050* (50 per Case)**

Macrobore Extension Set.
30", 5.0 mL priming volume. Male/Female
luer lock and pinch clamp.



▪ **Order #: ISSBC21* (50 per Case)**

Bifurcated Extension Set.
8", 1.5 mL priming volume. Male/Female
luer locks with back check valves.



▪ **Order #: ISB21* (50 per Case)**


Bifurcated Extension Set.
8", 2.1 mL priming volume. Male/Female
luer locks with pinch clamps.




▪ **Order #: ISM2* (50 per Case)**

Bifurcated Extension Set.
4", 0.3 mL priming volume. Male/Female
luer locks with slide clamps.


**Indicates special order, please contact representative for details*




- **Order #: ISMB21* (50 per Case)**
Bifurcated Minibore Extension Set.
5", 0.3 ml Male/Female Luer locks with pinch clamps.



- **Order #: ISMT31 (50 per Case)**
Minibore Trifurcated Extension Set.
8", 0.6 mL priming volume with Male/Female Luer locks and slide clamps.



- **Order #: B1025 (50 per Case)**
Minibore Trifurcated Extension Set.
2", 0.33 mL priming volume. Male/Female Luer locks with removable slide clamps.



- **Order #: CC-100 (100 per Case)**
Disposable Hydrogel Probe Cover

**Indicates special order, please contact representative for details*



NeoChild

(888) 887-NICU (6428)

www.neochild.com



NeoChild[®]